


## Introducing GLO Guided Light Optics Technology



### THE SCIENCE

Patented GLO technology combines optimal heat and light in a closed system mouthpiece to activate and accelerate the whitening gel formulation for safe, superior results with little to no sensitivity.

Healthy patients GLO!



### GLO DUAL WHITENING

**GLO professional whitening in the practice plus patient take-home kits**

Fast appointments. Superior results. High compliance. No lab fees. No trays. Little to no sensitivity. Clinically proven.

**For the best whitening experience – GLO!**



A three-time recipient of the prestigious Thomas Edison Award, the GLO Teeth Whitening Device, Everyday GLO Whitening Maintenance Pen, and GLO Solo G-Vials received silver, gold, and gold honors respectively for innovation in personal care.

Before GLO



After GLO



Before GLO



After GLO



Standard 32 minute treatment: **8 shades whiter**

Tetracycline Treatment: **16 shades whiter** in 4 weeks

Clinically Proven Whitening Results 8+ shades whiter – depending on starting shade.